



# Dipper™

## Urine Dipstick Control

- Simulates patient testing
- Liquid, made with real human urine
- 18 month unopened and 3 month (or 20 dips) open vial stability when stored at 2-8 C\*

### Analytes\*

Bilirubin	Creatinine	hCG	Leukocytes	Nitrite	Protein	Urobilinogen
Blood	Glucose	Ketones	Microalbumin	pH	Specific Gravity	

\*Refer to the package insert of currently available lots for specific analyte and stability claims

This is the one that started it all. We craft Dipper Urine Dipstick Control to monitor the performance of visual and instrumental readings of urine dipsticks by immersing the dipstick into the control – the same way you test your patient samples. We include values for all dipstick analytes plus microalbumin and creatinine; and qualitative results for hCG early pregnancy detection test methods. Dipper Urine Dipstick Control is designed for use with most urinalysis reagent strips. It can also be used for confirmatory tests and refractometry. Dip it, read it, done.

### ORDERING INFORMATION:

1440-01 2 Level Set..... 6 x 15 mL

1442-61 Level 2..... 6 x 15 mL



For More Information: **phone** +1.310.536.0006 **email** info@4qc.com **web** www.4qc.com