



QuanTtest Red

Total Protein Assay System

- For qualitative determination of protein in cerebrospinal fluid and urine for both manual and automated systems
- Reagents stable until the expiration date when stored properly*
- Standards and controls are stable until the expiration date when stored at 2-8 C*

**Refer to the package insert of currently available lots for specific analyte and stability claims*

We craft Pyrogallol Red Molybdate dye binding method to provide a simple, ready-to-use colorimetric method for total protein quantitation with greater linearity, using microliter samples of urine or cerebrospinal fluid (CSF) in manual or automated systems.

Assay kit consists of four reagent standards, reagent, and controls for both CSF and urine chemistry. Reagent standards include level 2 (25 mg/dL), 3 (50 mg/dL), 4 (100 mg/dL), and 5 (200 mg/dL).

ORDERING INFORMATION:

2210-02 Kit..... 200 Tests
5210-12 Reagent..... 1 x 250 mL



For More Information: **phone** +1.310.536.0006 **email** info@4qc.com **web** www.4qc.com